

EUROPEAN COMMISSION JOINT RESEARCH CENTRE

Directorate F - Health, Consumers and Reference Materials (Geel) Food and Feed Compliance



JRC F.5/CvH/MGH/AS/Ares

Evaluation Report on the Analytical Methods submitted in connection with the Application for Authorisation of a Feed Additive according to Regulation (EC) No 1831/2003

L-lysine monohydrochloride and L-lysine sulphate produced by fermentation with Corynebacterium glutamicum CGMCC 14498 (FAD-2020-0067; CRL/200043)



Evaluation Report on the Analytical Methods submitted in connection with the Application for Authorisation of a Feed Additive according to Regulation (EC) No 1831/2003

Dossier related to:	FAD-2020-0067 - CRL/200043
Name of Product:	<i>L-lysine monohydrochloride and L-lysine sulphate produced by fermentation with Corynebacterium glutamicum</i> CGMCC 14498
Active Agent (s):	L-lysine
Rapporteur Laboratory:	European Union Reference Laboratory for Feed Additives (EURL-FA) JRC Geel, Belgium
Report prepared by:	María José González de la Huebra
Report checked by: Date:	Stefano Bellorini 08/03/2021
Report approved by: Date:	Christoph von Holst 08/03/2021



EXECUTIVE SUMMARY

In the current application authorisation is sought under Article 4(1) for *L-lysine monohydrochloride* and *L-lysine sulphate* produced by fermentation with *Corynebacterium glutamicum* CGMCC 14498, under the category/functional group 3(c) 'nutritional additives'/'amino acids, their salts and analogues', according to Annex I of Regulation (EC) No 1831/2003. Authorisation is sought for all animal species.

According to the Applicant, *L-lysine monohydrochloride* has a minimum purity (mass fraction on dry matter basis) of 98.5 % and 78.8 % with respect to the *L-lysine* content, respectively, while *L-lysine sulphate* has a minimum purity of 73.0 % and 55.0 % with respect to the *L-lysine* content, respectively.

The two forms of the *feed additive* are intended to be added directly into *feedingstuffs* or through complementary feed, *premixtures* and *water*. However the Applicant did not propose any minimum or maximum content of *L-lysine* in *feedingstuffs*.

For the quantification of *lysine* in the *feed additive* the Applicant proposed the ring-trial validated method EN ISO 17180:2013 based on ion-exchange chromatography coupled to visible or fluorescence detection (IEC-VIS/FLD). This standard method does not distinguish between the salts of amino acids and cannot differentiate between enantiomers. It applies for products containing more than 10 % of amino acid. The following performance characteristics are reported: a relative standard deviation for repeatability (RSD_r) ranging from 0.7 to 1.7 % and a relative standard deviation for reproducibility (RSD_R) ranging from 1.5 to 2.5 %. In addition, the EURL identified the "L-lysine monohydrochloride monograph" of the Food Chemical Codex (FCC) for the identification of *L-lysine monohydrochloride* in the *feed additive* and the generic European Pharmacopoeia monograph (Ph. Eur. 20301) for the identification of sulphate ions in *feed additive*.

For the quantification of *lysine* in *premixtures, feedingstuffs* and *water* the Applicant suggested using the ring-trial validated European Union method (Commission Regulation (EC) No 152/2009) based on IEC coupled with photometric detection (IEC-VIS). This method, designed only for the analysis of amino acids in *premixtures* and *feedingstuffs*, does not distinguish between the salts and the amino acid enantiomers. The following performance characteristics were reported for the quantification of total *lysine*: RSD_r ranging from 2.1 to 2.8 % and RSD_R ranging from 3.0 to 6.7 %. Furthermore, as concluded in previous amino acids reports, the EURL also considers the IEC-VIS/FLD procedure described above as fit-for-purpose for the determination of *lysine* in *water*.

Based on the performance characteristics available, the EURL recommends for official control (i) the "L-lysine monohydrochloride monograph" of the Food Chemical Codex (FCC)



for the identification of *L-lysine monohydrochloride* in the *feed additive*; (ii) the European Pharmacopoeia monograph (Ph. Eur. 01/2008:20301) for the identification of sulphate ions in *feed additive*; (iii) the ring-trial validated method EN ISO 17180:2013 based on IEC-VIS/FLD to quantify free *lysine* in the *feed additive* and *premixtures* (containing more than 10 % *lysine*); (iv) the European Union method based on IEC-VIS for the quantification of *lysine* in *premixtures* and *feedingstuffs*; and (v) the ion-exchange chromatography methods coupled with post-column derivatisation and optical detection (IEC-VIS/FLD) or coupled with post-column derivatisation and photometric detection (IEC-VIS) for the quantification of *lysine* in *water*.

Further testing or validation of the methods to be performed through the consortium of National Reference Laboratories as specified by Article 10 (Commission Regulation (EC) No 378/2005), as last amended by Regulation (EU) 2015/1761) is not considered necessary.

KEYWORDS

L-lysine monohydrochloride and *L-lysine sulphate* produced by fermentation with *Corynebacterium glutamicum* CGMCC 14498, nutritional additives, amino acids, all animal species

1. BACKGROUND

In the current application authorisation is sought under Article 4(1) (authorisation of a new feed additive) for L-lysine monohydrochloride and L-lysine sulphate produced by Corynebacterium CGMCC fermentation with glutamicum 14498, under the category/functional group 3(c) 'nutritional additives'/'amino acids, their salts and analogues', according to Annex I of Regulation (EC) No 1831/2003. Authorisation is sought for all animal species [1]. These two forms of L-lysine produced by fermentation with different Corynebacterium glutamicum strains are currently authorised as feed additives under Commission Implementing Regulation (EU) 2019/1964 [2] and Commission Implementing Regulation (EU) 2020/1798 [3].

According to the Applicant, *L-lysine monohydrochloride* has a minimum purity (mass fraction on dry matter basis) of 98.5 % and 78.8 % with respect to the *L-lysine* content, respectively, while *L-lysine sulphate* has a minimum purity of 73.0 % and 55.0 % with respect to the *L-lysine* content, respectively. [4].

The *feed additives* are produced by fermentation with a genetically modified strain of *Corynebacterium glutamicum*. The production strain is deposited in the "Chinese General Microbiological Culture Collection Center" (CGMCC) under accession number CGMCC 14498.



The two forms of the *feed additive* are intended to be added directly into *feedingstuffs* or through complementary feed, *premixtures* and *water*. However in the case of *water* special care has to be taken to avoid any unbalance due to addition through *feed* and *water*. The Applicant did not propose any minimum or maximum content of *L-lysine* in *feedingstuffs* [5].

Note: The EURL has previously evaluated the analytical methods for the determination of *lysine* in the frame of several dossiers [6].

2. TERMS OF REFERENCE

In accordance with Article 5 of Regulation (EC) No 378/2005, as last amended by Regulation (EU) 2015/1761, on detailed rules for the implementation of Regulation (EC) No 1831/2003 of the European Parliament and of the Council as regards the duties and the tasks of the European Union Reference Laboratory concerning applications for authorisations of feed additives, the EURL is requested to submit a full evaluation report to the European Food Safety Authority for each application or group of applications. For this particular dossier, the methods of analysis submitted in connection with *L-lysine monohydrochloride* and *L-lysine sulphate* and their suitability to be used for official controls in the frame of the authorisation were evaluated.

3. EVALUATION

Description of the analytical methods for the determination of the active substance in the feed additive, premixtures, feedingstuffs and when appropriate water (section 2.6.1 of the dossier - Annex II of Commission Regulation (EC) No 429/2008)

For the quantification of *lysine* in the *feed additive* the Applicant proposed the ring-trial validated method EN ISO 17180:2013 - "Animal feeding stuffs - Determination of lysine, methionine and threonine in commercial amino acid products and premixtures" [7]. This standard method is based on the experimental protocol described in the EU method for the analysis of free amino acids (including *lysine*) [8]. It does not distinguish between the salts of amino acids and cannot differentiate between enantiomers. It applies for products containing more than 10 % of amino acid.

Free *lysine* is extracted with diluted hydrochloric acid and further diluted with a sodium citrate buffer. After addition of norleucine as internal standard, the amino acids are separated by High Performance Liquid Chromatography (HPLC) with an Ion Exchange Column (IEC). Free *lysine* is quantified either after post-column derivatisation with ninhydrine and visible (VIS) detection at 440 nm and 570 nm or by fluorescence detection (FLD) after post-column reaction with ortho-phthaldialdehyde with a detector excitation wavelength at 330 nm and



emission at 460 nm. The performance characteristics reported for the quantification of free *lysine* are listed in Table 1.

Based on the performance characteristics available, the EURL recommends for official control the EN ISO 17180:2013 method for the quantification of free *lysine* in the *feed additive* and *premixtures* (containing more than 10 % *lysine*).

For the quantification of *lysine* in *premixtures*, *feedingstuffs* and *water* the Applicant suggested using the European Union (EU) method dedicated for the determination of amino acids in feed [8]. This method was designed for the quantification of free (synthetic and natural) and of total (peptide-bound and free) amino acids in *premixtures* and *feedingstuffs*, using an amino acid analyser or IEC coupled with post-column derivatisation and VIS detection. It does not distinguish between the salts of amino acids and cannot differentiate between enantiomers.

Table 1: Method performance characteristics obtained in the frame of ring-trial validation studies				
(EN ISO 17180:2013 [7], European Union method [8] and EN ISO 13903:2005 [9]) for the				
determination of total lysine in the feed additive, premixtures and feedingstuffs.				

Ring-Trial	Matrix	<i>lysine</i> content g/kg	RSD _r %	RSD _R %
[7]	Feed Additive	459	0.8	2.3
	Premix 3	208	1.3	2.5
	Premix 4	168	1.3	2.3
	Premix 5	128	0.7	1.9
	Premix 6	123	1.7	2.1
	Premix 7	104	1.2	1.8
	Premix 8	102	1.2	1.5
	Premix 9	240	1.1	2.2
	Premix 10	233	0.8	1.8
	L-Lysine-HCl	760	0.9	1.8
[8]	Mixed pig feed	10	2.8	3.2
	Broiler compound	14	2.1	5.4
	Protein concentrate	48	2.4	3.0
	Premix	98	2.1	6.7
[9]	Poultry meal	3.6	3.1	9.9
	Broiler finisher feed	3.5	3.5	9.0
	Broiler starter feed	1.4	2.4	9.0
	Corn	0.3	3.1	13.1
	Fishmeal	4.2	2.8	7.9

RSD_n, RSD_R - relative standard deviation for repeatability and reproducibility, respectively



The free amino acids are extracted with diluted hydrochloric acid. Co-extracted nitrogenous macromolecules are precipitated with sulfosalicylic acid and removed by filtration. The solution is filtered and adjusted to pH 2.2. The amino acids are separated by IEC and determined by post-column derivatisation with ninhydrin and photometric detection at 570 nm. The procedure chosen for the determination of the total amino acids depends on the amino acids under investigation. *Lysine* can be determined in either oxidised or non-oxidised samples. Oxidation is performed at 0 °C with a performic acid/phenol mixture. The excess of oxidation reagent is decomposed with sodium disulfite. The oxidised or non-oxidised sample is hydrolysed with hydrochloric acid (6 mol/L) for 23 hours. The hydrolysate is adjusted to pH 2.2. The amino acids are separated by IEC and determined by post-column derivatisation with ninhydrin and photometric detection at 570 nm.

The EU method was ring-trial validated using four different matrices listed in Table 1. This method was further ring-trial validated by twenty-three laboratories, resulting in the EN ISO 13903:2005 method [9]. The performance characteristics reported for the quantification of total *lysine* are listed in Table 1. Furthermore, the following limits of quantification were reported for free *lysine* and total *lysine*: 0.04 and 0.3 g/kg *feedingstuffs*, respectively [9].

In the frame of stability studies the Applicant provided experimental data of the applicability of the EU method for the determination of *lysine* in *water* for both products [10].

Based on the performance characteristics available, the EURL recommends for official control the ring-trial validated EU method based on IEC-VIS to quantify *lysine* in *premixtures feedingstuffs* and *water*. Furthermore, as concluded in previous amino acids reports [6], the EURL also considers the IEC-VIS/FLD procedure described above as fit-for-purpose for the determination of *lysine* in *water*.

Methods of analysis for the determination of the residues of the additive in food (section 2.6.2 of the dossier - Annex II of Commission Regulation (EC) No 429/2008)

The evaluation of corresponding methods of analysis is not relevant for the present application.

Identification/Characterisation of the feed additive (section 2.6.3 of the dossier - Annex II of Commission Regulation (EC) No 429/2008)

The EURL found the Food Chemical Codex (FCC) "L-lysine monohydrochloride monograph", where different tests (including optical rotation) are used for the identification of *L-lysine monohydrochloride* in the *feed additive* [11] and thus recommends the Food Chemical Codex monograph for the identification of *L-lysine monohydrochloride*.

While for the identification of *L-lysine sulphate* in the *feed additive* a specific method is not available, the EURL found the generic European Pharmacopoeia monograph (Ph. Eur.



01/2008:20301) for the identification of the sulphate ions in *L-lysine sulphate* [12]. The EURL recommends for official control the European Pharmacopoeia monograph for the identification of the sulphate ions in the *feed additive*.

Further testing or validation of the methods to be performed through the consortium of National Reference Laboratories as specified by Article 10 (Commission Regulation (EC) No 378/2005, as last amended by Regulation (EU) 2015/1761), is not considered necessary.

4. CONCLUSIONS AND RECOMMENDATIONS

In the frame of this authorisation the EURL recommends for official control (i) the "L-lysine monohydrochloride monograph" of the Food Chemical Codex (FCC) for the identification of *L-lysine monohydrochloride* in the *feed additive*; (ii) the European Pharmacopoeia monograph (Ph. Eur. 01/2008:20301) for the identification of the sulphate ions in the *feed additive*; (iii) the ring-trial validated method EN ISO 17180:2013 based on ion-exchange chromatography coupled to visible or fluorescence detection (IEC-VIS/FLD) to quantify free *lysine* in the *feed additive* and *premixtures* (containing more than 10 % lysine); (iv) the European Union method based on IEC-VIS for the quantification of *lysine* in *premixtures* and *feedingstuffs*; and (v) the ion-exchange chromatography coupled with post-column derivatisation and photometric detection (IEC-VIS/FLD) or coupled with post-column derivatisation and photometric detection (IEC-VIS) methods for the quantification of *lysine* in *water*.

Recommended text for the register entry (analytical method)

For the identification of *L*-lysine monohydrochloride in the feed additive:

- Food Chemical Codex "L-lysine monohydrochloride monograph"

For the identification of sulphate in the *feed additive*:

- European Pharmacopoeia monograph 20301

For the quantification of *lysine* in the *feed additive and premixtures* (containing more than 10% *lysine*):

ion-exchange chromatography coupled with post-column derivatisation and optical detection (IEC-VIS/FLD) – EN ISO 17180

For the quantification of *lysine* in *premixtures* and *feedingstuffs*:

 ion-exchange chromatography coupled with post-column derivatisation and photometric detection (IEC-VIS), Commission Regulation (EC) No 152/2009 (Annex III, F)



For the quantification of *lysine* in *water*:

 ion-exchange chromatography coupled with post-column derivatisation and optical detection (IEC-VIS/FLD); or ion-exchange chromatography coupled with postcolumn derivatisation and photometric detection (IEC-VIS)

5. DOCUMENTATION AND SAMPLES PROVIDED TO EURL

In accordance with the requirements of Regulation (EC) No 1831/2003, reference samples of *L-lysine monohydrochloride* and *L-lysine sulphate* produced by fermentation with *Corynebacterium glutamicum* CGMCC 14498 have been sent to the European Union Reference Laboratory for Feed Additives. The dossier has been made available to the EURL by EFSA.

6. REFERENCES

- *Application, Reference SANTE/E5: FORW. APPL. 1831/0059-2020 & Annex I submission number 1600193863550-2670
- [2] Commission Implementing Regulation (EU) 2019/1964 of 26 November 2019 concerning the authorisation of L-lysine base, liquid, L-lysine monohydrochloride, liquid, L-lysine monohydrochloride, technically pure, and L-lysine sulphate as feed additives for all animal species O.J. L 307, 28.11.2019
- [3] Commission Implementing Regulation (EU) 2020/1798 of 30 November 2020 concerning the authorisation of L-lysine monohydrochloride produced by Corynebacterium glutamicum DSM 32932 and L-lysine sulphate produced by Corynebacterium glutamicum KFCC 11043 as feed additives for all animal species O.J. L 402, 01.12.2020
- [4] *Technical dossier, Section II: II.2.1.3. Qualitative and quantitative composition
- [5] *Technical dossier, Section II: II.5.1. Proposed mode of use in animal nutrition
- [6] EURL Evaluation Reports:

https://ec.europa.eu/jrc/sites/jrcsh/files/finrep-fad-2013-0045_l-lysine%20sulphate.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2016-0052_lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2018-0012_l-lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2018-0019_l-lysinehcl_sulphate.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2018-0028_l-lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2018-0037-lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2019-0015_l-lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2019-0015_l-lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2019-0014-lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2019-0014-lysine.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2019-00160028-lysinehcl.pdf https://ec.europa.eu/jrc/sites/jrcsh/files/finrep_fad-2019-00160028-lysinehcl.pdf

[7] EN ISO 17180:2013 - Animal feeding stuffs – Determination of lysine, methionine and threonine in commercial amino acid products and premixtures



- [8] Commission Regulation (EC) No 152/2009 of 27 January 2009 laying down the methods of sampling and analysis for the official control of feed, O.J. L 54, 26.02.2009
- [9] EN ISO 13903:2005- Animal feeding stuffs Determination of amino acids content
- [10] *Technical dossier, Section II: Annex 2.1.3a & Annex 2.1.3b
- [11] Food Chemical Codex monograph "L-lysine monohydrochloride", FCC 7 (2010), p.598
- [12] European Pharmacopoeia Monograph 01/2008:20301 Identification reactions of ions and functional groups – sulphates

*Refers to Dossier no: FAD-2020-0067

7. RAPPORTEUR LABORATORY & NATIONAL REFERENCE LABORATORIES

The Rapporteur Laboratory for this evaluation is the European Union Reference Laboratory for Feed Additives, JRC, Geel, Belgium. This report is in accordance with the opinion of the consortium of National Reference Laboratories as referred to in Article 6(2) of Commission Regulation (EC) No 378/2005, as last amended by Regulation (EU) 2015/1761.

8. ACKNOWLEDGEMENTS

The following National Reference Laboratories contributed to this report:

- Laboratori Agroalimentari, Departament d'Agricultura, Ramaderia, Pesca, Alimentació i Medi Natural. Generalitat de Catalunya, Cabrils (ES)
- Centro di referenza nazionale per la sorveglienza ed il controllo degli alimenti per gli animali (CReAA), Torino (IT)
- Państwowy Instytut Weterynaryjny, Pulawy (PL)
- Thüringer Landesanstalt für Landwirtschaft (TLL). Abteilung Untersuchungswesen. Jena (DE)
- Österreichische Agentur für Gesundheit und Ernährungssicherheit (AGES), Wien (AT)
- Ruokavirasto Helsinki¹ (FI)
- Instytut Zootechniki Państwowy Instytut Badawczy, Krajowe Laboratorium Pasz, Lublin (PL)
- Ústřední kontrolní a zkušební ústav zemědělský (ÚKZÚZ), Praha (CZ)
- Wageningen Food Safety Research² (WFSR) (NL)

¹ Name and address according to according COMMISSION IMPLEMENTING REGULATION (EU) 2015/1761: Elintarviketurvallisuusvirasto/Livsmedelssäkerhetsverket (Evira), Helsinki/Helsingfors .

² Name and address according to according COMMISSION IMPLEMENTING REGULATION (EU) 2015/1761: RIKILT Wageningen UR, Wageningen.