



JRC CRL - FCM Database / GC-MS

Sample Name: 2,4-bis(2,4-dimethylphenyl)-6-(2-hydroxy-4-octyloxyphenyl)-1,3,5-triazine
 Solvent: Hexane CAS: 002725-22-6
 Concentration: 100 µg/mL PM ref: 38885 [U.R.N. 132]



m/z	Abundance	Ion Intensity %	m/z	Abundance	Ion Intensity %
41.10	2165,0	6,56	377.10	4140,0	12,54
43.10	2555,0	7,74	378.10	2316,0	7,02
116.10	2329,0	7,06	397.30	3655,0	11,07
130.10	4604,0	13,95	508.40	2649,0	8,03
131.10	7450,0	22,57	509.40	33008,0	100,00
132.10	7542,0	22,85	510.40	11529,0	34,93
253.00	2998,0	9,08	511.40	2819,0	8,54
315.10	3220,0	9,76			

Spectrometer: HEWLETT PACKARD GC/MS 6890/5973
 Inlet system: capillary GC/MS
 Scan Range: 40-700 amu
 Source temperature: 230 °C

Flow: 1.2 mL/min
 Column: DB 17-HT (30 m x 0.25 mm x 0.15 µm)
 Programme temperature: 40 °C (3 min); 20 °C/min
 (350 °C); 350 °C (20 min)